

**When referring to this protocol, please cite:** Wang CC, Bajikar SS, Jamal L, Atkins KA, Janes KA. (2014) A time- and matrix-dependent TGFBR3–JUND–KRT5 regulatory circuit in single breast epithelial cells and basal-like premalignancies. *Nat Cell Biol*, 16, 345-56.

1. Dissociate  $\sim 10^6$  cells extensively with 250–500  $\mu$ l trypsin
  - The protocol is compatible with a variety of plating schemes and densities
  - *Fewer cells will not make it through the entire procedure*
  - *Complete trypsinization (as long as 30–45 min for some cells) is important to obtain a fully dissociated single-cell suspension*
2. Break up cell clumps by adding 1 ml PBS to each well and pipetting up and down several times.
3. Transfer cells from each well to a microcentrifuge tube and spin at 800 rcf for 3 min.
4. Resuspend cells in 1 ml PBS and spin at 800 rcf for 3 min.
5. For standard flow cytometry without live-dead discrimination, skip to Step #9.
  - *Live-dead discrimination is helpful for identifying non-apoptotic forms of cell death and for artifactual staining that could be observed in nonviable cells*
  - *The live-dead stain should not be used if a DAPI counterstain is planned or the FL1 channel is going to be used without compensation*
6. Resuspend cells in 1 ml PBS and add 1  $\mu$ l violet-fluorescent reactive dye (Invitrogen #L34955; dissolve each vial in 55  $\mu$ l DMSO and use within 1–2 weeks) to each sample.
  - *This dye can be used to exclude dead cells on flow cytometers with a violet laser line*
7. Incubate on ice for 30 min and spin at 800 rcf for 3 min.
8. Resuspend cells in 1 ml PBS and spin at 800 rcf for 3 min.
9. Resuspend cells in 160  $\mu$ l PBS.
  - *There is typically  $\sim 20$   $\mu$ l residual volume before adding the PBS, bringing the total volume to  $\sim 180$   $\mu$ l*
10. Spike in 20  $\mu$ l 37% PFA, vortex, and incubate at room temperature for 15 min.
11. Spin at 800 rcf for 3 min, and place samples on ice.
12. Aspirate the PFA solution and immediately permeabilize cells with 200  $\mu$ l ice-cold MeOH.
  - *Cells can be stored at  $-20^\circ\text{C}$  for days to weeks at this point.*
13. Change cells to a fresh tube
  - *Moving the cell suspensions to a fresh tube allows a better pellet to be formed at the first centrifugation, when a large fraction of cells can be lost*
14. Spin at 800 rcf for 3 min.
15. Wash cells with 500  $\mu$ l PBS + 0.1% Tween-20 (PBS-T), and spin at 800 rcf for 3 min.
16. Resuspend cells in 100  $\mu$ l 1 $\times$  Western blocking reagent (Roche #11921673001; aliquot 10 $\times$  stock and freeze at  $-20^\circ\text{C}$ ) in PBS-T and incubate for 30 min at room temperature.
  - *Western blocking reagent provides a more stringent block than goat serum or BSA and gives greater specificity to the staining; however, this block may be too stringent for some antibodies*
  - *An alternative low-stringency blocking solution is 1% BSA (use throughout all the subsequent staining steps)*
17. Spin at 800 rcf for 3 min.
18. Resuspend cells in 100  $\mu$ l of primary antibody solution diluted in 1 $\times$  Western blocking reagent in PBS-T and incubate for 1 hr at room temperature.
  - *Antibody dilutions must be determined empirically for each primary (1:50, 1:100, and 1:200 or 0.1 – 1  $\mu$ g is a good starting range of dilutions, but some will require up to ten-fold less antibody)*
19. Spin at 800 rcf for 3 min.
20. Wash cells with 500  $\mu$ l PBS-T, and spin at 800 rcf for 3 min.
21. Repeat Step #20.
22. Resuspend cells in 100  $\mu$ l of secondary antibody solution diluted in 1 $\times$  Western blocking reagent in PBS-T and incubate for 1 hr at room temperature.
  - *Antibody dilutions must be determined empirically for each secondary*

## General flow cytometry

Janes Lab Protocols

Entered by Kevin Janes

6/20/13

- *We have optimized the following secondaries:*
  - *Alexa 488-conjugated goat anti-mouse/rabbit (Invitrogen) for FL1 channel — 1:1000 dilution (0.2  $\mu$ g)*
  - *PE-conjugated goat anti-mouse/rabbit (Jackson ImmunoResearch) for FL2 channel — 1:200 dilution (0.25  $\mu$ g) PE is light sensitive and will fade over time with storage*
  - *PerCP-conjugated goat anti-mouse/rabbit (Jackson ImmunoResearch) for FL3 channel — 1:100 dilution (0.5  $\mu$ g) PerCP is very light sensitive... keep protected at all times!*
  - *Alexa 647-conjugated goat anti-mouse/rabbit (Invitrogen) for FL4 channel — 1:5000 dilution (0.04  $\mu$ g)*
- *PE is the brightest as for a single-channel experiment; Alexa 488 pairs with Alexa 647 with essentially zero compensation required, but prevents doublet discrimination on standard FACSCalibur instruments; Alexa 488 paired with PerCP allows doublet discrimination but will require slight compensation to account for bleedthrough*

23. Wash cells with 500  $\mu$ l PBS-T, and spin at 800 rcf for 3 min.

24. Repeat Step #23.

25. Resuspend cells in 200  $\mu$ l PBS-T and transfer to flow cytometry tubes (Falcon #352008).