

When referring to this protocol, please cite: Janes KA, Wang CC, Holmberg KJ, Cabral K, Brugge JS. (2010) Identifying single-cell molecular programs by stochastic profiling. *Nat Methods*, 7, 311-7.

1. Transport embedded samples on dry ice to the cryomicrotome
2. Place sample in the cryomicrotome box and equilibrate the box temperature to -24°C
 - *Higher box temperatures (e.g., -20°C) are acceptable but can cause distortion of the acini at the periphery of the Matrigel during sectioning*
 - *Slide rack must be inside the cryomicrotome box to store slides after sectioning*
 - *Keep the slides at room temperature*
3. Replace microtome blade (if disposable)
4. Spray microtome blade and anti-roll bar with a Kimwipe moistened with both ethanol and RNase Away
 - *RNase Away by itself will freeze to the blade*
 - *Be careful not to cut yourself when cleaning the blade*
5. Pop the embedded sample out of the cryomold, mount on the sectioning platform with Neg-50, and freeze until the sample has solidified on the sectioning platform
6. Mount the sectioning platform on the cryomicrotome and cut trim sections until ~ 3 mm of Matrigel is visible on the cutting surface
 - *Matrigel will appear as light pink against the white Neg-50; test slides can be prepared and checked under an inverted microscope to confirm that acini are on the slides*
7. Cut $8\ \mu\text{m}$ sections, adjusting the anti-roll bar to prevent the section from accordioning
 - *Do not graze sections (subsequent section will be $16\ \mu\text{m}$ instead of $8\ \mu\text{m}$)*
 - *Alternatively, one can omit the anti-roll bar and use the brush technique*
 - *$8\ \mu\text{m}$ sections are roughly one-cell thick to allow for the maximum recovery of biological material after microdissection*
8. Wick sections onto Superfrost plus slides (VWR #48311-703) and move the slide immediately to the slide rack inside the cryomicrotome box
 - *Do not dry the slide at room temperature (this will ruin the RNA integrity in the tissue)*
 - *To wick a second section per slide, warm the back of the slide with your finger for a few seconds before wicking (do not warm the part of the slide that contains the existing section)*
9. Continue cutting and wicking sections until good cuts are no longer possible
 - *One embedded sample typically yields 30-60 slides with two sections per slide*
 - *Do not return a half-sectioned sample to -80°C (the exposed sample will dehydrate and become damaged)*
10. Move the slide box containing the frozen sections to dry ice and store at -80°C for up to one month