

When referring to this protocol, please cite: Wang L, Brugge JS, Janes KA. (2011) Intersection of FOXO and RUNX1 gene-expression programs in single breast epithelial cells during morphogenesis and tumor progression. *Proc Natl Acad Sci*, 108, E803-12.

I. RNA isolation and purification

1. Aspirate the cell culture medium and lyse cells directly in 1 ml RNA STAT-60 (Tel-Test #CS-110) per $\sim 10^6$ cultured cells.
 - *This roughly corresponds to 1 ml per well for a 6-well plate or 0.5 ml for a 12-well plate (samples smaller than 12-well are not recommended)*
 - *Chamber slides containing 3D structures can be lysed in 0.5 ml RNA STAT-60 per chamber*
 - *Washing with PBS is not necessary and the delay can degrade unstable mRNAs*
 - *Note the expiration date on RNA STAT-60; the phenol will oxidize over time*
2. In a chemical fume hood, shear the denatured genomic DNA in each sample by pipetting the lysate up and down several times with a P200 micropipet.
 - *For 3D structures, use a P1000 to help dissolve the matrigel thoroughly*
3. After the lysate has become less viscous, transfer to a microcentrifuge tube and incubate at room temperature for 5 min to complete the dissociation of mRNAs from ribonuclear proteins.
 - *Lysates can be stored at this point for several weeks at -80°C*
 - *Make sure that the lysate does not touch the cap of the tube or else it can contaminate the later steps in the procedure; if the cap is contaminated, transfer to a fresh microcentrifuge tube*
4. Add 200 μl chloroform per ml of RNA STAT-60 and shake the samples vigorously for 15 sec.
 - *Many samples can be shaken in parallel by stacking two microcentrifuge racks on top of one another with the tubes in between*
 - *Chloroform causes the lysate to separate into aqueous and organic phases with denatured proteins at the interface*
 - *Chloroform degrades with light and with air over time, so small quantities should be purchased*
5. Centrifuge the samples for 15 min at max speed on a benchtop centrifuge in the cold room.
6. In a chemical fume hood, transfer the upper aqueous phase of the samples to a fresh microcentrifuge tube.
 - *$\sim 40\text{--}60\%$ of the volume of the RNA STAT-60 used for the lysis will be present in the aqueous phase; the volume of the aqueous phase depends upon the sample that was lysed*
 - *Be very careful not to pick up any of the organic phase or the protein interface at this step; do not touch anywhere on the microcentrifuge tube when extracting the very top of the aqueous phase*
7. Add 0.5 ml isopropanol per ml RNA STAT-60, mix well by vortexing, and incubate for 30 min at 4°C .
 - *Use molecular biology-grade isopropanol that is reasonably fresh (we have run into problems with yellow precipitates when "old" isopropanol is used here)*
8. Centrifuge the samples for 10 min at max speed on a benchtop centrifuge in the cold room.
 - *A white glassy pellet may be visible after this step, but for small starting numbers of cells, the pellet may be invisible*
 - *Orient the tubes in the centrifuge so that you can anticipate where the pellet will be formed*
9. Carefully aspirate the supernatant and wash the pellets with 1 ml 75% EtOH per ml RNA STAT-60
 - *It is not critical to remove all of the isopropanol and risk aspirating the pellet; simply aspirate so that less than 100 μl isopropanol remains*
10. Mix gently by inverting 5–6 times and centrifuge at 7500 rcf for 5 min at 4°C .
 - *Do not vortex vigorously or the pellet will break up and then not pellet well at the next step*
11. Carefully aspirate the supernatant and allow the pellet to air dry for 5-10 min at room temperature.
 - *Aspirate so that less than 100 μl EtOH remains and then pull off the remaining liquid with a micropipet and gel-loading tip; tip the microcentrifuge tube sideways and then draw off the liquid far away from the base of the tube where the pellet is located*
 - *If the pellet is drawn up by the micropipet, redispense the liquid into the tube and spin at max speed on a benchtop centrifuge for 15 sec to resediment the pellet*

- *Wet pellets can be warmed for a few seconds at 37°C to increase evaporation*
 - *Smaller or broken pellets will need to have the remaining liquid evaporated off at 42°C or, if 42°C is ineffective, 55°C*
12. Resuspend in 10 µl nuclease-free H₂O and incubate for 10 min at 37°C
 - *Resuspend each pellet as soon as the white precipitate has dried to a brownish hazy film at the base of the microcentrifuge tube*
 - *Do not overdry the pellet, because this will make it much more difficult to redissolve*
 - *When resuspending, make sure that the pipet tip does not have any residual liquid in it; if it does, that is a chunk of RNA that has not completely resuspended, and the sample must be pipetted up and down thoroughly to wash off this chunk*
 13. Quantify the RNA concentration on a Nanodrop and store at –80°C.
 - *After mixing the sample, pull 2 µl carefully off the top to avoid particulate material that will cause unstable absorbance readings*
 - *To confirm that the RNA has fully dissolved, measure each sample two separate times, removing the sample, pipetting it up and down in the tube extensively, and then re-adding it to the Nanodrop for the second reading—concentrations between separate measurements should be ±20% (if not, repeat until measurements stabilize)*
 - *The A260/A280 ratio will only be 1.6–1.7 because the sample was dissolved in H₂O rather than elution buffer*
 - *The A260/A230 ratio will be low because of residual guanidinium contamination from the RNA STAT-60; however, we have not seen this adversely affect the downstream first-strand synthesis for qRT-PCR*
 - *This procedure does not yield clean enough DNA for sequencing, microarrays, etc.*

II. TURBO DNase treatment

1. Remove residual genomic DNA with the TURBO DNA-free kit (Ambion # AM1907). Assemble the following reactions in PCR tubes:
 2, 1, or 0.4 µg total RNA (use a constant amount across all samples)
 2 µl 10× TURBO DNase buffer
 1 µl TURBO DNase
 Volume to 20 µl with nuclease-free H₂O
 - *Add as much RNA as you can based on the yield from the different samples that will be analyzed simultaneously*
2. Incubate for 30 min at 37°C on a PCR thermocycler (“DNase treatment” protocol on MyCycler)
3. Add 2 µl of DNase Inactivation Reagent and incubate for 5 min at room temperature, flicking the tubes 2-3 times during the inactivation to keep the Inactivation Reagent resuspended.
4. Spin in a PCR centrifuge for 2 min and then transfer 5 µl of the supernatant to a fresh PCR tube for first-strand cDNA synthesis
 - *This will correspond to 500, 250, or 100 ng DNase-treated RNA per first-strand reaction depending upon what starting amount was used in Step #1*
 - *Take an additional 2.5 µl of DNase-treated RNA from two samples you expect to be different; mix these together to provide a “no RT” control sample that checks for contaminating genomic DNA in RT-qPCR experiments*
 - *Also, put 5 µl of nuclease-free H₂O in a fresh PCR tube to serve as a “Blank” control sample for first-strand cDNA synthesis and RT-qPCR*

III. First-strand cDNA synthesis

1. To 5 µl DNase-treated RNA, add the following as a master mix:
 0.5 µl 50 µM oligo(dT)₂₄ primer
 0.5 µl 10 mM dNTPs
 0.5 µl nuclease-free H₂O
 6.5 µl total volume

RNA isolation-purification and first-strand cDNA synthesis

Entered by Kevin Janes

Janes Lab Protocols

12/28/16

- *We run the reverse transcription reaction at half the manufacturer's recommended volume, which is sufficient for 100 RT-qPCR experiments (see Janes_RTqPCR.pdf)*
2. Heat denature for 5 min at 65°C on a PCR thermocycler ("RNA denaturation" protocol on MyCycler) and then place on ice.
 3. Prepare the following as a master mix:
 - 2 µl 5× first-strand buffer
 - 0.5 µl 0.1 M DTT
 - 0.5 µl RNAsin Plus RNase inhibitor (Promega #N2615)
 - 0.5 µl Superscript III RT
 - 3.5 µl total volume
 - *Remember not to add this master mix to the no RT control. This sample must get a separate mix lacking Superscript III.*
 4. Mix thoroughly by pipetting and then incubate RT reactions at 50°C for 60 min, followed by heat denaturation at 70°C for 15 min on a PCR thermocycler ("SSIII rt" protocol on MyCycler).
 5. Transfer cDNA samples to a fresh microcentrifuge tube and store at -20°C