

Before starting:

- Have small-sample cDNA material as template for reamplification
- Have access to a 96-well plate or SPRIPlate ring magnet (Beckman #A32782 or equivalent)

I. Identification of optimum cycle threshold for reamplification

1. Set up a 19 μ l master mix-1 on ice
2. Add 1 μ l of 5 \times -diluted amplified cDNA to each well of an optically clear PCR strip tube (Bio-Rad #TLS0801 and #TCS0803) or hard-shell PCR plate (BioRad #HSP9601 and #MSB1001) (one well for each sample)
3. Add 19 μ l of the master mix to each well on ice and mix by pipetting
 - *Be sure to use barrier tips when mixing to avoid amplicon contamination of the pipette*
4. Run 40 cycles of qPCR:
 - i. 1 min at 94°C
 - ii. 2 min at 42°C
 - iii. 3 min at 72°C (fluorescence read at the end of the extension)
5. The optimum cycle for reamplification is when the first sample hits the middle of its exponential phase
 - *Going longer than the middle of the exponential phase will overamplify the sample and cause loss of quantitative accuracy*
 - *The Nextera XT kit requires only 1 ng of cDNA, meaning that most samples will have acceptable yield even if there is a great disparity in the apparent abundances of cDNA material for different samples*
 - *High-Fidelity gives the highest yield at the end of the reamplification (do not substitute AmpliTaq)*

II. Reamplification and purification of small-sample cDNA

1. Set up a 99 μ l master mix-2 on ice
2. Add 1 μ l of amplified cDNA to each well of a PCR strip tube or PCR plate (one well for each sample)
 - *Any thin-walled strip tube or plate will suffice at this step*
3. Add 99 μ l of the master mix to each well and mix by pipetting
 - *Be sure to use barrier tips when mixing to avoid amplicon contamination of the pipette*
4. Run the optimum number of cycles of PCR as determined in Part I:
 - i. 1 min at 94°C
 - ii. 2 min at 42°C
 - iii. 3 min at 72°C
5. Purify and concentrate the PCR reaction with AMPure XP beads (Beckman #A63880)
 - i. Equilibrate AMPure beads at room temperature for 30 min before starting
 - ii. Add 70 μ l AMPure beads per 100 μ l reamplification (0.7 \times volume), mix by pipetting 10 times, and incubate at room temperature for 15 min
 - *The lower volume of AMPure beads compared to the manufacturer's protocol enables selective purification of amplicon >300 bp in size and the removal of AL1-oligo(dT) concatemers*
 - iii. Magnetize the samples at room temperature for 5 min
 - *Depending on the magnet, some of the beads may stick to the upper walls of the tube; if using strip tubes, it is possible to tilt the strips and guide the beads to the bottom of the well*
 - *For the SPRIPlate ring magnet, the beads are supposed to stick to the walls of the tube*
 - iv. On the magnet, aspirate the supernatant with a gel-loading pipette tip, leaving 5 μ l in the well
 - *Leaving residual volume is important for plate magnets that pull the beads to the bottom of the well*
 - v. On the magnet, wash with 200 μ l of freshly prepared 80% EtOH for 30 sec and aspirate with a gel-loading pipette tip
 - *80% EtOH is defined here as 8 ml 100% EtOH + 2 ml H₂O (no voluming up)*
 - *Add the 80% EtOH gently to the beads to avoid disrupting them*

- vi. Repeat Step 3.v, making sure to remove all the residual EtOH
 - vii. Air dry the wells for 10 min at room temperature
 - *It is important to dry the beads thoroughly to improve recovery, as EtOH will reduce the solubility of the bead-bound amplicon*
 - *However, it is also important not to overdry the beads, as the amplicon will not elute off the bead surface; avoid a “cracked” appearance of dried beads on the tube walls*
 - viii. Resuspend the beads in 10 μ l of elution buffer (10 mM Tris-HCl, pH 8.5) and incubate for 2 min at room temperature
 - *When dissolving in this small of a volume, it is important to rinse all the beads to the bottom of the well for maximum recovery*
 - *Be sure to use barrier tips when mixing to avoid amplicon contamination of the pipette*
 - ix. Magnetize the samples at room temperature for 1 min
 - x. Transfer the supernatant to a fresh PCR strip tube or PCR plate
6. Repurify the PCR reaction with AMPure XP beads (Beckman #A63880)
- *This second round of purification improves the amplicon purity from 85% to >95%*
 - i. Equilibrate AMPure beads at room temperature for 30 min before starting
 - ii. Add 7 μ l AMPure beads per 10 μ l purified sample (0.7 \times volume), mix by pipetting 10 times, and incubate at room temperature for 15 min
 - *The lower volume of AMPure beads compared to the manufacturer’s protocol enables selective purification of amplicon >300 bp in size and the removal of AL1-oligo(dT) concatemers*
 - iii. Magnetize the samples at room temperature for 5 min
 - *Depending on the magnet, some of the beads may stick to the upper walls of the tube; if using strip tubes, it is possible to tilt the strips and guide the beads to the bottom of the well*
 - *For the SPRIPlate ring magnet, the beads are supposed to stick to the walls of the tube*
 - iv. On the magnet, aspirate the supernatant with a gel-loading pipette tip, leaving 5 μ l in the well
 - *Leaving residual volume is important for plate magnets that pull the beads to the bottom of the well*
 - v. On the magnet, wash with 200 μ l of freshly prepared 80% EtOH for 30 sec and aspirate with a gel-loading pipette tip
 - *80% EtOH is defined here as 8 ml 100% EtOH + 2 ml H₂O (no voluming up)*
 - *Add the 80% EtOH gently to the beads to avoid disrupting them*
 - vi. Repeat Step 6.v, making sure to remove all the residual EtOH
 - vii. Air dry the wells for 10 min at room temperature
 - *It is important to dry the beads thoroughly to improve recovery, as EtOH will reduce the solubility of the bead-bound amplicon*
 - *However, it is also important not to overdry the beads, as the amplicon will not elute off the bead surface; avoid a “cracked” appearance of dried beads on the tube walls*
 - viii. Resuspend the beads in 10 μ l of elution buffer (10 mM Tris-HCl, pH 8.5) and incubate for 2 min at room temperature
 - *When dissolving in this small of a volume, it is important to rinse all the beads to the bottom of the well for maximum recovery*
 - *Be sure to use barrier tips when mixing to avoid amplicon contamination of the pipette*
 - ix. Magnetize the samples at room temperature for 1 min
 - x. Transfer the supernatant to a fresh PCR strip tube or PCR plate
7. Use 1 μ l of the 2x purified amplicon for quantification by Qubit assay, and store the remainder at -20°C before Nextera XT tagmentation

III. Amplicon quantification by Qubit assay (Life Technologies #Q32850) on the CFX96 instrument

1. Prepare a seven-point standard curve by serially diluting the Qubit dsDNA BR Standard #2 (at 100 ng/ μ l) twofold in Qubit dsDNA Standard #1 (at 0 ng/ μ l) to yield 100, 50, 25, 12.5, 6.25, 3.1, and 1.6 ng/ μ l

Small-sample cDNA reamplification for RNA-seq

Janes Lab Protocols

Entered by Kevin Janes

07/01/15

- *Only 1 μ l of each serial dilution will be used in the assay*
- *Also be sure to include a blank of 100% Qubit dsDNA BR Standard #1*
2. Add 1 μ l of each serial dilution or 2x purified amplicon to the base of an optically clear PCR strip tube (Bio-Rad #TLS0801 and #TCS0803) or hard-shell PCR plate (BioRad #HSP9601 and #MSB1001) (one well for each sample)
3. Prepare Qubit Working Solution by diluting Qubit dsDNA BR Reagent 200-fold in Qubit dsDNA BR Buffer
 - *199 μ l of Buffer + 1 μ l of Reagent is enough for 8–9 samples*
4. Add 19 μ l of Qubit working solution to each well and mix by pipetting 10 times
5. Seal the strip tube or plate and read once on the CFX96 instrument at 22°C after a hold at 22°C for 2 min
 - *The Qubit detection reagent has an excitation maximum at 510 nm and an emission maximum at 527 nm, which is close enough to SYBR Green (Ex/Em = 497/520 nm) to use any qPCR detection instrument*
 - *Our saved protocol is named "QUBIT"*
6. Export the ZPCR file from the instrument and load into the CFX Manager Software
7. Once the PCRD file had been created, deselect all empty wells on the plate
 - *Having only the active wells selected makes data export easier*
8. Go to Settings > Baseline Setting > No Baseline Subtraction to turn off the default background subtraction settings
 - *It is crucial to export the raw RFU values or else the data will be uninterpretable*
9. Go to Export > Export All Data Sheets > Excel 2003
 - *The file ending in "End Point Results.xls" will contain the end RFUs from the instrument for the selected wells*
 - *Export > Custom Export will export only one file of RFUs if set up correctly, although all the wells of the 96-well plate will be exported*
10. Open the exported RFU data in Excel, perform linear regression on the standard curve, and back calculate the DNA concentrations of the 2x purified amplicon samples
11. Dilute a fraction of the 2x purified amplicon samples to 0.2 ng/ μ l in H₂O for tagmentation with the Nextera XT kit.

Master mix and buffer recipes

- **Master mix-1** Prepare fresh on ice
 - 2 μ l 10 \times High-Fidelity PCR buffer without Mg²⁺
 - 2.8 μ l 25 mM MgCl₂
 - 0.04 μ l 100 mM dATP
 - 0.04 μ l 100 mM dCTP
 - 0.04 μ l 100 mM dGTP
 - 0.04 μ l 100 mM dTTP
 - 0.1 μ l 20 mg/ml BSA
 - 0.06 μ l 15 mg/ml AL1 primer
 - 0.2 μ l High-Fidelity polymerase
 - 0.05 μ l 100 \times SYBR green
 - 13.63 μ l PCR-grade H₂O
 - 19 μ l total volume
- **Master mix-2** Prepare fresh on ice
 - 10 μ l 10 \times High-Fidelity PCR buffer without Mg²⁺
 - 14 μ l 25 mM MgCl₂
 - 0.2 μ l 100 mM dATP
 - 0.2 μ l 100 mM dCTP
 - 0.2 μ l 100 mM dGTP
 - 0.2 μ l 100 mM dTTP
 - 0.5 μ l 20 mg/ml BSA
 - 0.3 μ l 15 mg/ml AL1 primer
 - 1 μ l High-Fidelity polymerase
 - 72.4 μ l PCR-grade H₂O
 - 99 μ l total volume