

1. Plate cells at the desired density on 12 mm round No.1.5 coverslips (Thomas #1217N79) at the base of a 24-well plate.
  - *The DuoLink PLA Fluorescence protocol uses reactions of 40  $\mu$ l total volume on a surface of 1 cm<sup>2</sup>*
  - *Our protocol uses 25  $\mu$ l total volume per coverslip to preserve reagents*
  - *No. 1.5 coverslips provide the proper working distance for microscope objectives*
2. Aspirate the culture medium and wash cells with 500  $\mu$ l PBS.
3. Fix coverslips in freshly prepared 3.7% PFA and incubate for 15 min at room temperature.
  - *Stock solutions of 37% PFA should be prepared according to Wang et al. Meth. Enzymol. 85:514 (1982) and stored in single-use aliquots at -20°C*
  - *Fix slides in a fume hood and dispose of PFA as hazardous waste*
4. Wash coverslips 3  $\times$  5 min in PBS.
5. Block slides in a humid chamber for 1 hr at room temperature in blocking solution: 1 $\times$  Western Blocking Reagent (Roche #11921673001) diluted in PBS + 0.3% Triton X-100.
  - *Humid chambers can be constructed with an inverted pipette-tip box containing a paper towel saturated with water*
  - *Coverslips should be handled with fine forceps and lifted off the bottom of the 24-well plate with a curled hypodermic needle while the PBS is still in the well*
  - *25  $\mu$ l of blocking solution is sufficient for each coverslip*
  - *Place each coverslip face-down on a piece of Parafilm in the humid chamber*
6. Add primary antibody at the appropriate dilution in blocking solution overnight at room temperature.
  - *25  $\mu$ l of primary-antibody solution is sufficient for each coverslip*
  - *Antibody dilutions must be determined empirically, although 1:100 or 1:200 dilutions are common*
  - *Place each coverslip face-down on a piece of Parafilm*
7. Wash coverslips 2  $\times$  5 min in 1x Wash Buffer A at room temperature.
  - *Washes can be done in a clean 24-well plate*
8. Vortex PLUS and MINUS PLA probes and dilute each 1:5 in the DuoLink Antibody Diluent
  - *5  $\mu$ l PLUS probe + 5  $\mu$ l MINUS probe + 15  $\mu$ l Antibody Diluent per coverslip*
  - *Prepare as a master mix for all samples*
9. Add PLA probes and incubate in a humid chamber for 1 hr at 37°C.
10. Wash coverslips 2  $\times$  5 min in 1x Wash Buffer A at room temperature.
  - *Washes can be done in the same 24-well plate as before*
11. Dilute 5x DuoLink Ligation Buffer to 1x in ddH<sub>2</sub>O and dilute Ligase 1:40 in 1x Ligation Buffer.
  - *0.625  $\mu$ l Ligase in 25  $\mu$ l 1x Ligation Buffer per coverslip*
  - *Prepare as a master mix for all samples*
12. Add Ligase and incubate in a humid chamber for 30 min at 37°C.
13. Wash coverslips 2  $\times$  5 min in 1x Wash Buffer A at room temperature.
  - *Washes can be done in the same 24-well plate as before*
14. Dilute 5x Amplification Buffer to 1x in ddH<sub>2</sub>O and dilute Polymerase 1:80 in 1x Amplification Buffer.
  - *0.3125  $\mu$ l Polymerase in 25  $\mu$ l 1x Amplification Buffer per coverslip*
  - *Prepare as a master mix for all samples*
15. Add Polymerase and incubate in a humid chamber for 100 min at 37°C.
16. Wash coverslips 2  $\times$  10 min in 1x Wash Buffer B at room temperature.
  - *Washes can be done in the same 24-well plate as before but different wells on the plate*
17. Wash coverslips in 0.01x Wash Buffer B for 1 min.
  - *Dilute 1x Wash Buffer B to 0.01x in ddH<sub>2</sub>O immediately before use*
  - *Washes can be done in the same 24-well plate as before*
18. Mount coverslips with 10  $\mu$ l DuoLink InSitu Mounting Media with DAPI.
19. Seal the edge of the coverslip with nail polish and allow to air dry.
20. Remove residual salts with a Kimwipe and air dry.
21. Store slides at 4°C until imaging.

## Proximity ligation assay on coverslips

Janes Lab Protocols

Entered by Kevin Janes

7/18/22

### Buffer recipes

- 1x Wash Buffer A (for 1 L total volume):

One Wash Buffer A pouch

1000 ml ddH<sub>2</sub>O

Store a 4°C long term and bring working aliquot to room temperature before use

- 1x Wash Buffer B (for 1 L total volume):

One Wash Buffer B pouch

1000 ml ddH<sub>2</sub>O

Store a 4°C long term, bring working aliquot to room temperature before use, and dilute to 0.01x at room temperature immediately before use