

When referring to this protocol, please cite: Wang L, Brugge JS, Janes KA. (2011) Intersection of FOXO and RUNX1 gene-expression programs in single breast epithelial cells during morphogenesis and tumor progression. *Proc Natl Acad Sci*, 108, E803-12.

1. Plate 293T cells at 400,000 cells/well on 6-well plates for 24 hr.
 - *Cells should be ~30% confluent on the day of transfection*
 - *Do not package virus with 293T cells that have been cultured for more than two months since thawing from liquid N₂; we have found that packaging efficiency drops dramatically after extended culturing*
2. One hour before transfection, aspirate medium and refeed cells with 2 ml growth medium (DMEM + 10% FBS + pen/strep)
 - *Do not wash cells with PBS, as this will cause most of the cells to detach*
 - *Refeeding before transfection increases transfection efficiency by ~2-fold*
3. For amphotropic retroviruses, prepare the following calcium phosphate solution for each well of the 6-well plate:
 - 1.25 µg retroviral vector (pBabe, pLNCX, etc.)
 - 1.25 µg pCL amphi (retroviral packaging vector)
 - 10 µl 2.5 M CaCl₂
 - Volume to 100 µl with 0.1× TE (pH 7.6)
 - *This “double transfection” protocol prepares much higher titer retroviruses than GPG, Phoenix, etc. cells that stably express the packaging genes*
 - *We generally transfect two well per retroviral construct to have enough virus for triple infection if needed (above solution can be scaled up without problems)*

For amphotropic lentiviruses, prepare the following calcium phosphate solution:

- 1.25 µg lentiviral vector (pLKO.1, etc.)
- 0.75 µg psPAX2 (lentiviral gag-pol packaging vector)
- 0.5 µg pMD.2G (VSV-G envelope protein)
- 10 µl 2.5 M CaCl₂
- Volume to 100 µl with 0.1× TE (pH 7.6)
- *Using the psPAX2-pMD.2G packaging vectors appears to improve virus titers by ~2-fold compared to D8.2-pCMV-VSVG packaging vectors*

For lentiviral reporters, prepare the following calcium phosphate solution:

- 0.25 µg pTRF.1 lentiviral reporter vector
- 1.75 µg pFIV-34N (lentiviral gag-pol packaging vector)
- 0.5 µg pVSV-G envelope vector
- 10 µl 2.5 M CaCl₂
- Volume to 100 µl with 0.1× TE (pH 7.6)
- *pFIV-34N and pVSV-G were cloned from the pPACKF1 plasmid mix from Systems Biosciences #LV100A-1*

4. To each 100 µl of 2× DNA-CaCl₂ mixture, add 100 µl of 2× HEPES-buffered saline at room temperature, and mix by pipetting
5. Let the precipitate stand for 1 min, add the entire 200 µl dropwise to the well of a 6-well plate, and gently agitate the plate to mix the precipitate with the medium
 - *Prepare the precipitates individually, because transfection efficiency is reported to drop quickly after 1 min*
 - *Removing the medium and adding the precipitate directly to the 293T cells does not increase transfection efficiency in my hands*
 - *Calcium phosphate precipitates will appear as small black speckles on the inverted microscope; it is good to confirm their presence to ensure that the precipitation was effective*
6. Incubate the precipitates with the cells for 4–6 hrs at 37°C, then aspirate the medium and carefully replace with 1 ml of growth medium

Amphotropic virus preparation in 293T cells by calcium phosphate

Entered by Kevin Janes

Janes Lab Protocols

2/19/09

- *The purpose of this medium replacement is to refeed the 293T cells, not remove the precipitates*
- 7. Collect supernatant at 48 hr for retroviruses and 24+48 hr for lentiviruses
- 8. Pass virus solutions through a 0.45 μm filter and store at 4°C for short-term storage (days) or -80°C for long-term storage
 - *For low m.o.i. infections, we have found that freshly prepared virus gives the highest titers*

Buffer recipes

- **2.5 M CaCl₂**
36.76 g CaCl₂
Volume to 100 ml in H₂O
Sterilize by autoclaving and store at room temperature
- **1x TE (pH 7.6)**
158 mg Tris-Cl
29 mg EDTA
Volume to 90 ml and pH to 7.6 with NaOH
Volume to 100 ml in H₂O and dilute to 0.1x in H₂O
Sterilize by autoclaving
Store at room temperature
- **2x HEPES-buffered saline**
800 mg NaCl
27 mg Na₂HPO₄•2H₂O
1.2 g HEPES
Volume to 90 ml in H₂O and pH to 7.05 with NaOH (0.5 N)
Volume to 100 ml in H₂O.
Sterilize by passing through a 0.45 μm filter
 - *HEPES solutions cannot be autoclaved*Store 1 ml and 5 ml aliquots at -20°C